WASHING MILK CONTAINERS

1. RINSE

Use <u>lukewarm</u> water. DO NOT rinse with hot water. Rinse off dirt and milk residue.

2.WASH

Use <u>hot</u> water. Add liquid detergent and bleach or a dry chlorinated detergent. Brush all surfaces. Scrub off remaining milk residue. Keep water above 49° (120°F) at all times.

3.RINSE

Use <u>warm</u> water. Add acid. Rinse containers. Do not rinse off the acid solution. Leave it on the bottles and pails while they dry. [May use an acid-sanitizer, also.]

4.DRY

Allow the bottles and pails to drain and dry. Do not stack pails inside each other. Do not sit pails upside down on a concrete floor.

> Sam Leadley, Calf & Heifer Management Specialist Shirley Macmillan, United Kingdom Editor <u>smleadley@yahoo.com</u> <u>www.atticacows.com</u>
> © Attica Vet. Assoc. 2017 All Rights Reserved